

# SYNGUARD®

## NITRILE EXAM GLOVES

### FEATURES

- Excellent softness and wearing fitness
- Finger textured for a reliable, consistent grip
- Beaded cuff ensures easy donning and prevents roll down
- Compliant with food handling requirements
- Latex-free, powder-free and odor-free



PRODUCT INFORMATION			
Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Available Sizes	XS,S,M,L,XL	Color	Blue, White, Violet, Black
Glove Length	≥240mm	Freedom from Holes (Inspection level I)	AQL 1.5
Powder Content	Powder-Free	Sterile	No
Finger Thickness(mm)	≥0.05	Palm Thickness (mm)	≥0.05
Force at break	≥6N before and after aging	Tensile strength	≥14MPA before and after aging
Protection against Bacteria and Fungi	Pass	Protection against 40% Sodium Hydroxide (K) 30% Hydrogen peroxide (P) 37% Formaldehyde (T)	>30min
Protection against Viruses	Pass		
Dimensions(±10mm)	XS/70;S/80;M/95; L/110;XL/120	Surface	Micro-textured on fingers

ORDERING INFORMATION						
Size		XS	S	M	L	XL
SKU	Blue	SNBE10013	SNBE10014	SNBE10015	SNBE10016	SNBE10017
	White	SNBE10023	SNBE10024	SNBE10025	SNBE10026	SNBE10027
	Violet	SNBE10033	SNBE10034	SNBE10035	SNBE10036	SNBE10037
	Black	SNBE10043	SNBE10044	SNBE10045	SNBE10046	SNBE10047

PACKAGING AND STORAGE		
Packaging	100pcs/box,10boxes/ctn	200pcs/box,10boxes/ctn
Box size	230*125*60mm	245*130*80mm
Carton Dimension	315*258*245mm	415*270*258mm
Shelf Life	5 Years (recommended)	
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.	

### STANDARDS AND CERTIFICATIONS

Compliant with EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009  
 Compliant with EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016, EN ISO 21420